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

Bronnen:

University of Reading: investigating photosynthesis

For investigation 1: Light and dark

- Hydrogen carbonate indicator solution (enough for each student / group to fill 2 bijou bottles, plus spare)
- Bijou bottles and lids (2 per student / group)
- 1 bright lamp (e.g. 42W portable CFL or 150W halogen floodlamp with heat sink) or several smaller desk lamps fitted with 100W filament or 100W equivalent energy saving (CFL) bulbs (see notes on lighting)
- Black sugar paper cut into strips to wrap around the bijou bottles

For investigation 2: Light intensity (distance from lamp)

- Hydrogen carbonate indicator solution (enough for each student / group to fill 2 bijou bottles, plus spare)
- Bijou bottles and lids (2 per student / group)
- 1 bright lamp (e.g. 42W portable CFL or 150W halogen floodlamp with heat sink) or several smaller desk lamps fitted with 100W filament or 100W equivalent energy saving (CFL) bulbs (see notes on lighting)
- Black sugar paper cut into strips to wrap around the bijou bottles
- 50cm rulers

For investigation 3: Light colour/wavelength

- Hydrogen carbonate indicator solution (enough for each student / group to fill 5 bijou bottles)
- Bijou bottles and lids (5 per student / group)
- 1 bright lamp (e.g. 42W portable CFL or 150W halogen floodlamp with heat sink) between 2-6 students or several smaller desk lamps fitted with 100W filament or 100W equivalent energy saving (CFL) bulbs (see notes on lighting)
- green, red, blue and clear acetate films, cut to wrap around bijou bottles (one of each per student / group)
- Black sugar paper cut into strips to wrap around the bijou bottles (one per student / group)

For accurate filters, we recommend Primary Red, Primary Green and Bray Blue. These can be obtained from Lee Filters (see suppliers section).

For investigation 4: Number of algal balls or concentration of algae in each ball

- Hydrogen carbonate indicator solution (enough for students to fill 5 bijou bottles each)
- Bijou bottles and lids (5 per student)
- 1 bright lamp (e.g. 42W portable CFL or 150W halogen floodlamp with heat sink) between 2-6 students or several smaller desk lamps fitted with 100W filament or 100W equivalent energy saving (CFL) bulbs (see notes on lighting)