

Developing a HyGenChip for detection of the genotoxic potency of hydrophilic contaminants



Antoine Karengera

antoine.karengera@wetsus.nl

## **Motivation**

Knowledge gaps related to hydrophilic chemicals in groundwater and surface water are a critical challenge in the management of water quality [1,2]. Currently, there is no existing proper method to concentrate and chemically analyze low concentrations of water soluble pollutants. Most of these contaminants remain unseen, unmeasured, and largely unknown. Despite their low concentrations, these substances present a potential risk to the drinking water chain and environmental water quality (Figure1). This especially is the case for very sensitive endpoints such as genotoxicity and endocrine disruption. This project aims at development of a bioassay to quantify the genotoxic potency of hydrophilic compounds. By exposing very small organisms to water with hydrophilic compounds, the molecular response in the organisms can be used as an indication of the genotoxic potency of compound mixtures in the water (Figure 2).

## **Technological challenge**

The challenge is to design an efficient and cheap test that does not require specialists to apply and interpret. The bioassay, HyGenChip, employs very small invertebrate organisms as sentinel species. Hydrophilic compounds, invisible by chemical analyses, will leave their signature in those animals. The focus is on genotoxic contaminants that damage DNA and induce metabolic and DNA repair responses. The most indicative gene expressions for genotoxic effects will be validated with standard single compounds and mixtures. With selected gene markers a MAGPIX assay (Figure 3 & 4) will be developed. The sample preparation needed for the MAGPIX assay is far less time consuming than for classical gene expresion analyses. The assay can then be performed and interpreted in normally equipped laboratories.



Figure 2. Graphical Abstract of the project



Figure 3. Schematic representation summerizing the concept of a DNA microarray method. This assay provides the possibility to monitor the expression levels of thousands of genes simultaneously [3].



Figure 4. Schematic representation of multiplex assay with the Magpix system for gene expression profing [4].

## **Research goals**

1. Develop an easily applicable method and approach 'HyGenChip'



*Figure 1. Schematic representation of some possible sources of hydrophilic contaminants and their pathways to groundwater and surface water*  for detection and quantification of the genotoxic potency of hydrophilic pollutants;

- To assess the genotoxic potency of hydrophilic contaminants in surface water (environmental relevance) and groundwater (source of raw drinking water);
- 3. Develop a simple and easy procedure to apply and interpret the approach for genotoxicological assessment.
- [1] Loos, R., et al. Water Res 47 (2013): 6475-6487
- [2] Lopez, B. et al. Sci. Total Environ. 518–519 (2015): 562–573
- [3] Nuwaysir, E.F. Et al. Mol Calcinog 24 (1999): 153-159
- [4] S. Angeloni, R. Cordes, S. Dunbar et al. (2016), xMAP Cookbook, 3<sup>rd</sup> edition, 2014, http://www.luminexcorp.com



This project has received funding from the European Union's Horizon 2020 research and innovation programme under the Marie Skłodowska-Curie grant agreement No 665874

## www.wetsus.eu www.wageningenur.nl/mae

Antoine Karengera MSc,. prof.dr. Tinka Murk, dr. Inez Dinkla, ing. Marije IJszenga

